

Highlights:

- **High Sensitivity Protocol**
Can detect as little as 0.17% Herculex™ I corn or WideStrike™ Cotton in a 2 hour assay
- **Rapid Protocol**
Screens single seed or leaf samples for the presence of Cry1F protein in a 1 hour assay

Contents of Kit:

- 1 antibody-coated plate
- Cry1F Positive Control
- Cry1F Enzyme Conjugate
- 1 packet of Buffer Salts
- Substrate
- Stop Solution



Prepare wash buffer and extraction solutions

Catalog Number AP 016

Intended Use

The EnviroLogix QualiPlate™ Kit for Cry1F is designed for the qualitative laboratory detection of Cry1F endotoxin in corn and cotton, seed, leaf and bulk seed/grain samples. Two assay protocols are presented: The High Sensitivity Protocol will detect the Cry1F protein found in 0.17% Herculex I corn or 0.17% WideStrike cotton (1 positive seed in a 600 seed sample) and requires two hours to run. The Rapid Protocol (one hour total) is intended for use in screening individual seed or leaf punch samples for the presence of Cry1F.

How the Test Works

The QualiPlate Kit is a “sandwich” Enzyme-Linked ImmunoSorbent Assay (ELISA).

In the test, sample extracts are added to test wells coated with antibodies raised against Cry1F protein. Any Cry1F present in the sample extract binds to the antibodies, and is then detected by addition of enzyme (horseradish peroxidase)-labeled Cry1F antibody.

After a simple wash step, the results of the assay are visualized with a color development step; color development is proportional to Cry1F concentration in the sample extract.

Lighter color = Lower concentration

Darker color = Higher concentration

Items Not Provided

- distilled or deionized water for preparing Wash and Extraction Buffer
- Tween-20 (Sigma P-1379, or equivalent) for preparation of Extraction Buffer
- glass bottles or flasks plus graduated cylinder with 1 liter capacity for preparation and storage of Wash and Extraction Buffer
- Waring laboratory blender (model 31BL91 or equivalent), glass jar adapter (Eberbach # E8495) and appropriate size glass Mason jars for ground seed samples
- Snap-cap tubes and pestles for extraction of leaf samples (EnviroLogix Cat# ACC 002, 100/package), optional
- centrifuge capable of 5000 x g, optional
- disposable tip, adjustable air-displacement pipettes which will measure 50 and 100 microliters (µL), preferably of multi-channel configuration
- marking pen (indelible)
- tape or Parafilm®
- timer
- microtiter plate reader or strip reader
- wash bottle, or microtiter plate or strip washer
- racked dilution tubes for loading samples into the plate with a multi-channel pipette (optional)
- orbital plate shaker (optional)

USDA Websites

www.gipsa.usda.gov/fgis/handbook/gihbk1_inspec.aspx

USDA Grain Inspection Handbook, Book 1, Grain Sampling. This document provides a comprehensive overview of recommended sampling guidelines for static lots and grain streams. It reviews the various types of equipment and strategies that can be used to obtain a representative grain sample from different types of containers.

www.gipsa.usda.gov/fgis/biotech/sample2.htm

Guidance document entitled Sampling for the Detection of Biotech Grains, which provides important statistical sampling considerations when testing for the presence of biotech grains. It covers the basis for making probability determinations in accepting lots based upon different assumptions with respect to sample size, number of samples, sample preparation, etc.

www.gipsa.usda.gov/fgis/biotech/sample1.htm

Practical Application of Sampling for the Detection of Biotech Grains. This one-page application guide provides a table that gives sample sizes for selected lot concentrations and probability of rejecting the specified concentrations. It also provides a formula for making the calculation for other combinations.

www.gipsa.usda.gov/fgis/biotech/samplingplan1.xls

This website provides a simple to use Sample Planner (29k Excel Spreadsheet). The planner allows you to enter different assumptions in terms of sample size, number of samples, acceptable quality level and to determine the probability of accepting lots with given concentration levels. It also plots the probabilities in graph form for easy interpretation. Specific data can be saved for documentation and future analyses.

Preparation of Solutions

Wash Buffer:

Add the contents of the packet of **Buffer Salts** (phosphate buffered saline, pH 7.4 - Tween 20) to 1 liter of distilled or deionized water and stir to dissolve. Store refrigerated when not in use; warm to room temperature prior to assay. If more wash buffer is needed, order item # P-3563 from Sigma Chemical Co. (St. Louis, MO), or prepare the equivalent.

Extraction Buffer:

Add 0.5 mL Tween-20 to 100 mL of prepared Wash Buffer, and stir to dissolve. Store refrigerated when not in use; warm to room temperature prior to assay.

Cry1F Positive Control:

The Positive Control is used as provided in the Rapid Protocol, but must be diluted 1:4 for use in the High Sensitivity Protocol. Prepare this dilution just prior to running the assay: mix 50 μ L Cry1F Positive Control with 150 μ L Extraction Buffer for each set of duplicate wells to be filled.

Sample Preparation

Note: It is recommended that the user prepare known negative and positive seed or leaf samples to be run in every assay as controls, in addition to the kit Positive Control.

Ground Grain/Seed:

This protocol requires that a small sample (20 to 50 grams) be analyzed. It is essential that this sample be well mixed and representative of the larger bulk. The test will detect 0.17% Herculex I corn or WideStrike cotton (one positive seed in a sample of 600 seeds). Bulk cottonseed samples should be delinted.

NOTE: Thorough mixing of the bulk grain sample and determination of an appropriate sampling plan are critical to the results of this testing, and are the responsibility of the user of this test kit. The USDA/GIPSA has prepared several guidance documents to address the issues involved in obtaining representative grain samples from static lots—such as trucks, barges, and railcars—and for taking samples from grain streams.

Sampling plans should be chosen that best meet the needs of both the buyer and seller in terms of acceptable risks. Increasing the number of kernels in the sample and taking multiple samples will increase the likelihood of obtaining representative samples, and maximize the probability of detecting any contamination in the grain lot. For further information on USDA/GIPSA guidelines for obtaining representative samples and assessing detection probabilities for biotech grain, see the websites listed to the left.

It is the responsibility of the user to ensure proper sampling and thorough mixing prior to analysis. Once representative samples have been obtained from the truck or container, they can be reduced in size using a splitter and uniformly ground and mixed. **The finer the grind, the faster and more efficient the extraction.**

1. For 600 seed samples, grind in a “Mason” jar (32 oz. for corn, 16 oz. for delinted cottonseed) on a blender at high speed for 1 minute. Shake jar to mix, then repeat the grinding a second time. Thoroughly clean the grinding equipment between samples to prevent cross-contamination.
2. Weigh at least 20 grams of ground sample into a jar or cup.




Extract grain sample



Centrifuge to clarify grain extract



Leaf punch using microtube cap



Add Extraction Buffer Blank, Negative and Positive Control, and each sample extract to the plate

3. Add 50 mL of Extraction Buffer to each 20 gram corn sample, or 60 mL to each cotton sample. For all other grain sample sizes, add Extraction Buffer at the rate of 2.5 mL per gram of corn or 3.0 mL per gram of cotton. Cap and shake vigorously by hand or vortex for 20-30 seconds. Let stand at room temperature for 1 hour to extract. Mix again at the end of the hour.
4. Clarify the extracts by centrifuging at 5000 x g for 5 minutes. Alternatively, allow them to settle out for at least 10 minutes. Insert a pipette tip below any floating lipid layer and above the pellet to remove the clarified sample. Dispensing particles into the test plate can cause false positive results.
5. Use the High Sensitivity Protocol to test these sample extracts.

Single Seeds or Leaf Punches:

Individual seeds:

1. Crush seeds: Seeds may be crushed by any number of methods, from hammers or pliers in a bag or tube, to 48-well seed crushers, to bead-beater type grinders. Whatever the method used, take extreme care not to cross-contaminate between seed samples.
2. Add 1 mL of Extraction Buffer to each crushed seed. Mix for at least 30 seconds. For best results, allow to extract for an hour, mixing again at the end of that time. If seeds are thoroughly crushed, this extraction time can be reduced. Allow extracts to settle completely. Dispensing particles into the test plate can cause false positive results.


Leaf testing:

1. Take a single leaf punch of approximately 5 to 10 millimeters diameter, using a paper punch or micro-tube cap. Mash the leaf tissue with a pestle matched to the micro-tube, or disrupt via another method. The extraction efficiency of whatever method used will vary proportionately with the amount of tissue disruption performed.
2. Add 0.25 to 0.5 mL of Extraction Buffer per leaf punch. Mix for at least 30 seconds, and allow particles to settle. Take extreme care not to cross-contaminate between leaf samples. Dispensing particles into the test plate can cause false positive results.

Use the Rapid Protocol to test seed and leaf extracts.

How to Run the Assay


- Read all of these instructions before running the kit.
- Allow all reagents to reach room temperature before beginning (at least 30 minutes with un-boxed plates and reagents at room temperature - do not remove plates from bag with desiccant until they have warmed up).
- Organize all Controls, clarified sample extracts, and pipettes so that Step 1 can be performed in 15 minutes or less. The use of a multi-channel pipette is strongly recommended for all reagent additions.
- Use the well identification markings on the plate frame as a guide when adding the samples and reagents. In a qualitative assay, the Blank (BL), Positive Control (PC) in duplicate wells, and 92 sample extracts (S) in single wells may be run on one plate. (See the Qualitative Assay Example Plate Layout - Figure 1A).
- Choose the assay protocol that fits the testing needs.



Mix plate and incubate



Add Enzyme Conjugate



Mix plate and incubate



Bottle Wash method



Strip Plate Wash option

HIGH SENSITIVITY PROTOCOL

The **High Sensitivity Protocol** will detect 0.17% Herculex I corn or WideStrike cotton in ground grain/seed, and requires two hours of total assay incubation time. Dilute the Cry1F Positive Control 1:4 in Extraction Buffer for this protocol.

1. Add **50 μ L** of **Extraction Buffer Blank (BL)**, **50 μ L** of **diluted Positive Control (PC)**, and **50 μ L** of each **sample and user-prepared control extract (S)** to their respective wells, as shown in the Example Plate Layout (Figure 1A).

NOTE: In order to minimize setup time it is strongly recommended that a multi-channel pipette be used in steps 1, 5, 8, and 10.

2. Thoroughly mix the contents of the wells by moving the plate in a rapid circular motion on the benchtop for a full 20-30 seconds. Be careful not to spill the contents!
3. Cover the wells with tape or Parafilm to prevent evaporation and incubate at ambient temperature for **30 minutes**. If an orbital plate shaker is available, shake plate at 200 rpm.
4. After incubation, carefully remove the covering and vigorously shake the contents of the wells into a sink or other suitable container. Flood the wells completely with **Wash Buffer**, then shake to empty. Repeat this wash step three times. Alternatively, perform these four washes (300 μ L/well) with a microtiter plate or strip washer. Slap the inverted plate on a paper towel to remove as much liquid as possible.
5. Add **50 μ L Cry1F-enzyme Conjugate** to each well. Thoroughly mix the contents of the wells, as in step 2.
6. Cover the wells with tape or Parafilm to prevent evaporation and incubate at ambient temperature for **1 hour**. If an orbital plate shaker is available, shake plate at 200 rpm.
7. Wash the plate as described in step 4.
8. Add **100 μ L of Substrate** to each well.
9. Thoroughly mix the contents of the wells, as in step 2. Cover the wells with new tape or Parafilm and **incubate for 30 minutes at ambient temperature**. Use orbital shaker if available.

Caution: Stop Solution is 1.0N Hydrochloric acid. Handle carefully.

10. Add **100 μ L of Stop Solution** to each well and mix thoroughly. This will turn the well contents yellow.

NOTE: Read the plate within 30 minutes of the addition of Stop Solution.

RAPID PROTOCOL

The **Rapid Protocol** is used to screen single seed or leaf samples for presence of Cry1F and requires one hour of total assay incubation time. Use the Cry1F Positive Control without dilution for this protocol.

1. Add **50 μ L Cry1F-enzyme Conjugate** to each well. Immediately add **50 μ L** of **Extraction Buffer Blank (BL)**, **50 μ L** of **Positive Control (PC)**, and **50 μ L** of each **sample extract and user-prepared control extract (S)** to their respective wells, as shown in the Example Plate Layout (Figure 1A).



Slap inverted plate on towel to remove as much liquid as possible



Complete protocol and add Stop Solution



Read plates in a Plate Reader within 30 minutes of the addition of Stop Solution

NOTE: In order to minimize setup time it is strongly recommended that a multi-channel pipette be used in steps 1, 5 and 7.

2. Thoroughly mix the contents of the wells by moving the plate in a rapid circular motion on the benchtop for a full 20-30 seconds. Be careful not to spill the contents!
3. Cover the wells with tape or Parafilm to prevent evaporation and incubate at ambient temperature for **45 minutes**. If an orbital plate shaker is available, shake plate at 200 rpm.
4. After incubation, carefully remove the covering and vigorously shake the contents of the wells into a sink or other suitable container. Flood the wells completely with **Wash Buffer**, then shake to empty. Repeat this wash step three times. Alternatively, perform these four washes (300 μL /well) with a microtiter plate or strip washer. Slap the inverted plate on a paper towel to remove as much liquid as possible.
5. Add **100 μL of Substrate** to each well.
6. Thoroughly mix the contents of the wells, as in step 2. Cover the wells with new tape or Parafilm and **incubate for 15 minutes at ambient temperature**. Use orbital shaker if available.

Caution: Stop Solution is 1.0N Hydrochloric acid. Handle carefully.

7. Add **100 μL of Stop Solution** to each well and mix thoroughly. This will turn the well contents yellow.

NOTE: Read the plate within 30 minutes of the addition of Stop Solution.

How to Interpret the Results

Spectrophotometric Measurement

1. Set the wavelength of the microtiter plate reader to 450 nanometers (nm). (If it has dual wavelength capability, use 600, 630 or 650 nm as the reference wavelength.)
2. Set the plate reader to **blank** on the **Extraction Buffer Blank** wells (this should automatically subtract the mean optical density (OD) of the Blank wells from each control and sample OD). If the reader cannot do this, it must be done manually.

General test criteria:

- The mean OD of the BLANK wells should not exceed 0.15.
- The mean, blank-subtracted OD of the Positive Control wells should be at least 0.15.
- The coefficient of variance (%CV) between the duplicate Positive Control wells should not exceed 15%:

$$\%CV = \frac{\text{std. deviation of OD's}}{\text{mean Pos.Ctl. OD}} \times 100$$

If the results of an assay fail to meet these criteria, consult EnviroLogix' Technical Service for suggestions on improving the test when you repeat the assay.

Calculate the Positive Control Ratio

Divide the OD of each sample extract by the mean OD of the Positive Control wells. This number is the "Positive Control Ratio".

Interpret the Qualitative Results

Ground corn or cottonseed samples

If the Positive Control Ratio calculated for a sample is less than 0.5, the ground corn contains less than 0.17% Herculex I corn or WideStrike cotton.

If the Positive Control Ratio of a sample is greater than or equal to 0.5, the sample contains 0.17% or greater Herculex I corn or WideStrike cotton.

NOTE: This test is to be used qualitatively only, with yes/no results at 0.17% positive. For information on testing at different cutoff levels, please contact EnviroLogix' Technical Service.

Single leaf and seed samples:

If the Positive Control Ratio calculated for a sample is less than 1.0, the sample does not contain Cry1F at the levels normally found in Herculex I corn or WideStrike cotton.

If the Positive Control Ratio of a sample is greater than or equal to 1.0, the sample is Herculex I corn or WideStrike cotton.

Leaf and seed samples are by their nature either 100% positive or 100% negative. Any low level positive results from single seed or leaf samples must be due to either some form of sample cross-contamination (stray particles or dust, leaf residue on leaf punch, etc.) or can be caused by transfer of particulate matter from leaf or seed extracts into the assay wells. If there is any question of the latter occurring, re-extraction and re-testing is recommended.

Figure 1A. Example of a typical Qualitative assay setup.

	1	2	3	4	5	6	7	8	9	10	11	12
A	BL	S7	S15	S23	S31	S39	S47	S55	S63	S71	S79	S87
B	PC	S8	S16	S24	S32	S40	S48	S56	S64	S72	S80	S88
C	S1	S9	S17	S25	S33	S41	S49	S57	S65	S73	S81	S89
D	S2	S10	S18	S26	S34	S42	S50	S58	S66	S74	S82	S90
E	S3	S11	S19	S27	S35	S43	S51	S59	S67	S75	S83	S91
F	S4	S12	S20	S28	S36	S44	S52	S60	S68	S76	S84	S92
G	S5	S13	S21	S29	S37	S45	S53	S61	S69	S77	S85	BL
H	S6	S14	S22	S30	S38	S46	S54	S62	S70	S78	S86	PC

Precautions and Notes

- Store all QualiPlate Kit components at 4°C to 8°C (39°F to 46°F) when not in use.
- Do not expose QualiPlate Kit components to temperatures greater than 37°C (99°F) or less than 2°C (36°F).
- Allow all reagents to reach ambient temperature (18°C to 27°C or 64°F to 81°F) before use.
- Do not use kit components after the expiration date.
- Do not use reagents or test well strips from one Plate Kit with reagents or test well strips from a different Plate Kit.

- Do not expose **Substrate** to **sunlight** during pipetting or while incubating in the test wells.
- The assay has been optimized for use with the protocol provided in the kit. Deviation from this protocol may invalidate the results of the test.
- The results generated through the proper use of this kit reflect the condition of the working sample directly tested. Extrapolation as to the condition of the originating lot from which the working sample was derived should be based on sound sampling procedures and statistical calculations which address random sampling effects, non-random seed lot sampling effects, and assay system uncertainty. A negative result obtained when properly testing the working sample does not necessarily mean the originating lot is entirely negative for the analyte or protein in question.
- As with all tests, it is recommended that results be confirmed by an alternate method when necessary.
- Observe any applicable regulations when disposing of samples and kit reagents.
- Use caution to prevent sample-to-sample cross-contamination with samples, fluids, or disposables.



**For Technical Support
Contact Us At:**

EnviroLogix
500 Riverside Industrial
Parkway
Portland, ME 04103-1486
USA

Tel: (207) 797-0300
Toll Free: 866-408-4597
Fax: (207) 797-7533

e-mail:
info@envirologix.com

website:
www.envirologix.com



LIMITED WARRANTY

EnviroLogix Inc. (“EnviroLogix”) warrants the products sold hereunder (“the Products”) against defects in materials and workmanship when used in accordance with the applicable instructions for a period not to extend beyond a product’s printed expiration date. If the Products do not conform to this Limited Warranty and the customer notifies EnviroLogix in writing of such defects during the warranty period, including an offer by the customer to return the Products to EnviroLogix for evaluation, EnviroLogix will repair or replace, at its option, any product or part thereof that proves defective in materials or workmanship within the warranty period.

ENVIROLOGIX MAKES NO OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. The warranty provided herein and the data, specifications and descriptions of EnviroLogix products appearing in EnviroLogix published catalogues and product literature are EnviroLogix’ sole representations concerning the Products and warranty. No other statements or representations, written or oral, by EnviroLogix’ employees, agents or representatives, except written statements signed by a duly authorized officer of EnviroLogix Inc., are authorized; they should not be relied upon by the customer and are not a part of the contract of sale or of this warranty.

EnviroLogix does not warrant against damages or defects arising in shipping or handling, or out of accident or improper or abnormal use of the Products; against defects in products or components not manufactured by EnviroLogix, or against damages resulting from such non-EnviroLogix made products or components. EnviroLogix passes on to customer the warranty it received (if any) from the maker thereof of such non-EnviroLogix made products or components. This warranty also does not apply to Products to which changes or modifications have been made or attempted by persons other than pursuant to written authorization by EnviroLogix.

THIS WARRANTY IS EXCLUSIVE. The sole and exclusive obligation of EnviroLogix shall be to repair or replace the defective Products in the manner and for the period provided above. EnviroLogix shall not have any other obligation with respect to the Products or any part thereof, whether based on contract, tort, strict liability or otherwise. Under no circumstances, whether based on this Limited Warranty or otherwise, shall EnviroLogix be liable for incidental, special, or consequential damages.

This Limited Warranty states the entire obligation of EnviroLogix with respect to the Products. If any part of this Limited Warranty is determined to be void or illegal, the remainder shall remain in full force and effect.

This test kit has been validated and approved by Dow AgroSciences for detection of the Cry1F protein expressed in Herculex I corn and WideStrike cotton.

*Parafilm is a registered trademark of American Can Corporation
Herculex and WideStrike are trademarks of Dow AgroSciences LLC
Tween is a registered trademark of Uniqema, a business unit of ICI Americas Inc.
EnviroLogix, the EnviroLogix logo, and QualiPlate are trademarks of EnviroLogix Inc.*

© EnviroLogix 2015



Safety Data Sheet
According to OSHA 29CFR 1910.1200

SECTION 1. Identification of the substance/mixture and of the company/undertaking

1.1 Product identifier
Trade name: **Wash Buffer Salts**
Part number: 50-0091, 10099

1.2 Relevant identified uses of the substance or mixture and uses advised against application of the substance or the preparation: Laboratory chemicals

1.3 Details of the supplier of the safety data sheet
Manufacturer/Supplier: EnviroLogix Inc., 500 Riverside Industrial Pkwy, Portland ME 04103, USA (207) 797-0300 (207) 797-0300 Technical Service

1.4 Emergency telephone number:

SECTION 2. Hazards identification

2.1 Classification of the Substance or Mixture:
Classification according to OSHA 29CFR 1910.1200 (Hazard Communication): Not a hazardous substance or mixture

2.2 Label Elements:
None required according to 29CFR 1910.1200

Other indications: None

2.3 Additional information:
No other information

SECTION 3. Composition/information on ingredients

3.1 Mixture: Potassium acid phosphate
Synonyms: PBS

Hazardous Components	Chemical name	CAS No	EC No	Amount (%)	Classification
	Potassium Chloride	7447-40-7	231-211-8	1-5 %	Aquatic Acute 3; Aquatic Chronic 3, H412

Based on the amount of hazardous ingredients in this product, it is not considered hazardous according to 29CFR 1910.1200

SECTION 4. First aid measures

4.1 Description of first aid measures:
After inhalation: Supply fresh air, consult doctor in case of breathing difficulties.
After skin contact: Flush skin with plenty of water for at least 15 minutes. Remove contaminated clothing. Seek medical attention if irritation develops.
After eye contact: Rinse opened eye for several minutes under running water. Seek medical attention if irritation develops.
After swallowing: If swallowed, consult with medical staff or poison control center to determine if any immediate response or follow up actions are recommended. Never give anything by mouth to an unconscious person.

4.2 Most important symptoms and effects, both acute and delayed:
None

4.3 Indication of any immediate medical attention and special treatment needed:
No special treatment is required

SECTION 5. Firefighting measures

5.1 Extinguishing media:
Suitable extinguishing agents: CO₂, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.

5.2 Special hazards arising from the substance or mixture:
Carbon oxides, Oxides of Phosphorus, Potassium, Sodium, Hydrogen Chloride gas

5.2 Advice for firefighters:
Wear protective equipment appropriate for fire conditions including respiratory protective gear

SECTION 6. Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures:
Use PPE, avoid dust formation, ensure adequate ventilation, avoid breathing dust

6.2 Environmental precautions:
Prevent further leakage or spillage if safe to do so. Do not let product enter drains. Discharge to the environment must be avoided.

6.3 Methods and material for containment and clean up:
Pick up and arrange disposal without creating dust. Sweep up and shovel. Keep in suitable closed containers for disposal

6.4 Reference to other sections:
For safe handling refer to Section 7; For information on PPE refer to Section 8. For disposal, refer to Section 13.

SECTION 7. Handling and storage

7.1 Precautions for safe handling:
Practice good chemical hygiene when handling. Avoid contact with eyes, skin and clothing. Prevent formation of dust.

7.2 Conditions for safe storage, including any incompatibilities:
Keep containers closed, store in a dry, well ventilated space.

7.3 Specific end use(s):
Apart from the uses mentioned in section 1.2, no other end uses are stipulated.

SECTION 8. Exposure controls/personal protection

8.1 Control parameters:
Components with workplace control Parameters: Contains no substances with occupational exposure limit values

8.2 Exposure controls:
8.2.1 Appropriate engineering controls: Ensure eyewash and safety shower are nearby; provide ventilation if necessary

8.2.2 Personal Protective Equipment:
Eyes: Safety glasses with side shields, goggles. Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166 (EU). Eye and face protection equipment are described by OSHA (US) in 29CFR1910.133. Do not wear contact lenses when working with chemicals

Hands: Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands. The selected protective gloves have to satisfy the specifications of EU Directive 89/686/EEC and the standard EN 374 derived from it.

Respiratory protection: Appropriate respiratory protection should be determined according to local conditions using risk analysis protocols. An approved disposable air purifying particulate respirator may be used as a backup to engineering controls. Always use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

Body: Use body protection relative to its type and amount of material being handled

8.2.3 Environmental controls: Sweep or wipe up spills, do not allow into sewers or drains

SECTION 9. Physical and chemical properties

9.1 Information on basic physical and chemical properties:

a) Appearance:	White powder.
b) Odor:	None
c) Odor Threshold:	No data available
d) pH:	7-8
e) Melting point/freezing point:	No data available
f) Boiling point/boiling range:	No data available
g) Flash point:	No data available
h) Evaporation rate:	No data available
i) Flammability (solid, gaseous):	No data available
j) Upper/lower flammability or explosive limits:	No data available
k) Vapor pressure:	No data available
l) Vapor density:	No data available
m) Relative density:	No data available
n) Solubility(ies):	Water soluble
o) Partition coefficient: n-Octanol/water:	No data available
p) Auto-ignition temperature:	No data available
q) Decomposition temperature:	No data available
r) Viscosity:	No data available
s) Explosive properties:	No data available
t) Oxidizing properties:	No data available

9.2 Other information: No further relevant information available.

SECTION 10. Stability and reactivity

10.1 Reactivity: No data available

10.2 Chemical stability: Stable under normal recommended storage conditions.

10.3 Possibility of hazardous reactions: No data available

10.4 Conditions to avoid: No data available

10.5 Incompatible materials: Strong oxidizing agents and strong acids.

10.6 Hazardous decomposition products: No data available

SECTION 11. Toxicological information

Acute toxicity: No data available

Inhalation: No data available

Dermal: No data available

Skin corrosion/irritation: No data available

Serious eye damage: No data available

Respiratory or skin sensitization: No data available

Mutagenicity and toxicity for reproduction: No data available

Carcinogenicity: No component of this product at levels greater than 0.1 % is identified as probable, possible, or confirmed human carcinogen by IARC, ACGIH, NTP, or OSHA.

SECTION 12. Ecological information

12.1 Toxicity: No data available

12.2 Persistence and degradability: No data available

12.3 Bio accumulative potential: No data available

12.4 Mobility in soil: No data available

12.5 Results of PBT and vPvB assessment: Not available as a chemical safety assessment, not required/not conducted.

12.6 Other adverse effects: No data available

SECTION 13. Disposal considerations

Dispose of excess or unused product in accordance with Local, State and Federal regulations. Contact a licensed professional waste disposal service to dispose of this material.

SECTION 14. Transport information

14.1 UN-number (DOT, ADR, ADN, IMDG, IATA): Not dangerous goods

14.2 UN proper shipping name (DOT, ADR, ADN, IMDG, IATA): Not dangerous goods

14.3 Transport hazard classes (DOT, ADR, ADN, IMDG, IATA): Not applicable

14.4 Packing group (DOT, ADR, IMDG, IATA): Not applicable

14.5 Environmental hazards: Not applicable

14.6 Special precautions for user: Not applicable

14.7 Transport in bulk according to Annex II of MARPOL 73/78: Not applicable

SECTION 15. Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

US Federal Regulations
SARA Section 302 (Extremely Hazardous Substances): Not listed
Clean Air Act: Not listed
Clean Water Act: Not listed
OSHA: Not listed

US State Regulations
Massachusetts Right to Know: Disodium Hydrogenorthophosphate CAS No 7558-79-4 Rev Date: 2007-03-01
California Prop. 65 Components: Contains no chemicals known to cause cancer, birth defects, or reproductive harm.

15.2 Chemical Safety Assessment: Not carried out

SECTION 16. Other information

Hazard Code
H412 Harmful to aquatic life with long lasting effects

This information is true based on our present knowledge. However, EnviroLogix makes no representation of its accuracy or completeness. Persons receiving this information must exercise their independent judgment in determining the product's safety and suitability for its intended use. This document shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

EFES Department
EnviroLogix Inc.



Material Safety Data Sheet
OSHA 29CFR 1910.1200

SECTION 1. Identification of the substance/mixture and of the company/undertaking

1.1 Product identifier	Stop Solution
Trade name:	L.O.N HCl
Synonyms:	10825, 10827, 10828, 11193, 11776 (XGDD007)
Part number:	Laboratory chemicals
1.2 Relevant identified uses of the substance or mixture and uses advised against application of the substance / the preparation :	
1.3 Details of the supplier of the safety data sheet	Envirol ogix Inc., 500 Riverside Industrial Pkwy, Portland ME, 04103, USA Phone: (207) 7974300
1.4 Emergency telephone number:	(207) 797-0300 Technical Service

SECTION 2. Hazards identification

2.1 Classification of the substance or mixture	Hazard Classes
Classification according to OSHA 29 CFR 1910.1200	Metal Corrosive (Cat. 1) H290 Skin Irritation (Cat 2) H315 Serious Eye damage (Cat. 1) H318
2.2 Label elements	
Labeling according to OSHA 29CFR 1910.1200	
Hazard pictograms :	
Signal word :	Warning
Hazard statements:	H290 May be corrosive to metals H315 Causes skin irritation H318 Causes serious eye damage
Precautionary statements:	P281 Use personal protective equipment as required P302 + P352 IF ON SKIN: Wash with plenty of soap and water. P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do. Continue rinsing.
2.3 Other Statements	None

6.3 Methods and material for containment and cleanup:	Absorb in paper towel and discard in appropriate waste. Clean with water afterwards. Large spills may be neutralized with dilute solutions of sodium carbonate or calcium oxide.
6.4 References to other sections:	For safe handling refer to Section 7. For information on PPE refer to Section 8. For disposal refer to Section 13.

SECTION 7. Handling and storage

7.1 Precautions for safe handling:	Practice good chemical hygiene when handling. Avoid contact with eyes, skin, and clothing.
7.2 Conditions for safe storage, including any incompatibilities:	Store in tightly closed, non-metal container, in a corrosive compatible area. Prevent direct sunlight and heat. Store in well aired storage rooms.
7.3 Specific end use(s):	Apart from the uses mentioned in section 1.2., no other specific uses are stipulated.

SECTION 8. Exposure controls/personal protection

8.1 Exposure limits:	Components with limit values that require monitoring at the workplace:									
	<table border="1"> <thead> <tr> <th>Hydrogen Chloride</th> <th>European (Commission directive 96/94)</th> <th>USA (OSHA)</th> </tr> </thead> <tbody> <tr> <td></td> <td>8hr TWA = 5 ppm (7.5 mg/m³)</td> <td>Ceiling Limit = 5 ppm (7.5 mg/m³)</td> </tr> <tr> <td></td> <td>STEL = 10 ppm (15 mg/m³)</td> <td></td> </tr> </tbody> </table>	Hydrogen Chloride	European (Commission directive 96/94)	USA (OSHA)		8hr TWA = 5 ppm (7.5 mg/m ³)	Ceiling Limit = 5 ppm (7.5 mg/m ³)		STEL = 10 ppm (15 mg/m ³)	
Hydrogen Chloride	European (Commission directive 96/94)	USA (OSHA)								
	8hr TWA = 5 ppm (7.5 mg/m ³)	Ceiling Limit = 5 ppm (7.5 mg/m ³)								
	STEL = 10 ppm (15 mg/m ³)									
8.2 Exposure Controls:	Facilities using this mixture should be equipped with an eyewash and safety shower. Use general or local exhaust ventilation to keep airborne concentrations below permissible exposure limits.									
8.2.1 Engineering controls										
8.2.2 General protective and hygienic measures:	The usual precautionary measures should be adhered to when handling chemicals.									
Eye Protection:	Safety glasses with side shields, goggles. Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166 (EU). Eye and face protection regulations are described by OSHA (US) in 29CFR1910.133. Do not wear contact lenses when working with chemicals.									
Hand Protection:	Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands. The selected protective gloves have to satisfy the specifications of EU Directive 89/686/EEC and the standard EN 374 derived from it.									
Breathing Equipment:	Appropriate respiratory protection should be determined according to local conditions using risk analysis protocols. An approved disposable air purifying particulate respirator may be used as a backup to engineering controls. Always use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).									
8.2.3 Environmental exposure controls:	Contain spills, do not allow into environment									

SECTION 3. Composition/information on ingredients

3.2 Mixture	Aqueous solution 1N Hydrochloric Acid (1N HCl, 3% HCl)			
Chemical name	Amount (%)	CAS No		Classification According to OSHA 29CFR 1910.1200
		EC No		
Hydrochloric acid	1-4 %	7647-01-0		Hazard Classification May be Corrosive to Metals H290 Causes Skin Irritation H315 Causes Serious Eye Damage H318
		231-595-7		

SECTION 4. First aid measures

4.1 Description of first aid measures	
After inhalation :	In case of inhalation: Remove to fresh air. If not breathing give artificial respiration. Get medical attention immediately.
After skin contact :	In case of skin contact: Remove contaminated clothing and shoes immediately. Wash affected area with mild soap or detergent for at least 10 minutes or until no evidence of chemical remains.
After eye contact :	In case of eye contact, immediately flush eyes with plenty of water for at least 15 minutes. Lifting eyelids occasionally, until no evidence of chemical remains. Get medical attention immediately.
After swallowing :	In case of ingestion, DO NOT induce vomiting unless directed to do so by medical personnel. Never give anything by mouth to an unconscious person. Call a physician immediately.
4.2 Most important symptoms and effects, both acute and delayed:	May cause skin irritation and eye damage
4.3 Indication of any immediate medical attention and special treatment needed:	DO NOT use sodium bicarbonate in an attempt to neutralize the acid.

SECTION 5. Firefighting measures

5.1 Extinguishing media:	CO ₂ , extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
5.2 Special hazards arising from the substance or mixture:	Hydrogen Chloride gas
5.3 Advice for firefighters:	Wear protective gear appropriate for fire conditions including respiratory protective gear.

SECTION 6. Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures:	In the case of spilled mixture wear gloves to prevent skin contact. In the case of a large spill, additional protection is recommended.
6.2 Environmental precautions:	Do not discharge mixture to sewer system or waterways.

SECTION 9. Physical and chemical properties

9.1 Information on basic physical and chemical properties:	
a) Appearance:	Clear liquid, colorless to slight yellow.
b) Odor:	Pungent (slight)
c) Color Threshold:	No Data Available
d) pH:	pH 1
e) Melting point/freezing point:	No Data Available
f) Boiling point/Boiling range:	No Data Available.
g) Flash point:	Not applicable.
h) Evaporation rate:	0.36 (Water) compared with n-Butyl Acetate = 1
i) Flammability (solid, gaseous):	No Data Available
j) Upper/lower flammability or explosive limits:	No Data Available
k) Vapor pressure:	No Data Available
l) Vapor density:	No Data Available
m) Relative density:	No Data Available
n) Solubility(ies):	Fully miscible, water.
o) Partition Coefficient: n-Octanol/water:	No Data Available
p) Auto-ignition temperature:	No Data Available
q) Decomposition temperature:	No Data Available
r) Viscosity:	No Data Available but should be similar to that of water
s) Explosive properties:	No Data Available.
t) Oxidizing properties:	No Data Available
9.2 Other information:	No further relevant information available.

SECTION 10. Stability and reactivity

10.1 Reactivity:	No data available
10.2 Chemical Stability:	Stable under normal temperatures and pressures.
10.3 Possibility of hazardous reactions:	Under normal conditions of storage and use, hazardous reactions will not occur.
10.4 Conditions to avoid:	No specific data
10.5 Incompatible materials:	Metals, Alkali metals, bases, Amines.
10.6 Hazardous decomposition products:	Under normal conditions of storage and use, hazardous decomposition products should not be produced.

SECTION 11. Toxicological information

Information on Toxicological Effects			
Acute effects (toxicity tests):	7647-01-0 HCl	Effect Dose	Species
	Acute oral toxicity	LD50=900mg/kg	rabbit
	Acute dermal toxicity	No data	
	Acute inhalative toxicity	LC50 = 3124 mg/L	rat
Sensitization:	No sensitizing effects known		
CMR (carcinogenicity, mutagenicity and toxicity for reproduction) effects:	No CMR effects.		
Additional toxicological information:	No Additional Information		

SECTION 12. Ecological information

12.1 Toxicity:				
	Aquatic toxicity (1N HCl)	Effect dose	Exposure time	Species
	Acute fish toxicity	LC50=826 mg/L	96h	Lepomis idus
	Acute daphnia toxicity	No data		
	Acute algae toxicity	No data		