

Highlights:

- Test Cry34Ab1 in corn seed lot purity in 1 hour

Contents of Kit:

- 10 Cry34Ab1 antibody-coated solid plates
- Positive Control
- Cry34Ab1 Antibody-Enzyme Conjugate
- Substrate

Materials Not Provided

- PBS/0.05% Tween-20 Wash Buffer, pH 7.4
- distilled or deionized water for preparing Wash and Extraction Buffers
- Tween® 20 (Sigma cat# P1379), or equivalent for Extraction Buffer
- Stop Solution
- user-supplied controls: seed and/or leaf extracts from known negative and positive samples (optional)
- EnviroLogix Tissue Extraction Kit (ACC 002) or other suitable equipment for taking and extracting leaf punch samples
- equipment for pulverizing seeds or leaves
- disposable tip, adjustable air-displacement multi-channel pipettes which will measure 50 and 100 microliters (µL)
- marking pen (indelible)
- tape or Parafilm®
- timer
- microtiter plate reader with 450 nm filter
- wash bottle, or microtiter plate or strip washer
- orbital plate shaker (optional)

Catalog Number AP 054 NW V10

Intended Use

The EnviroLogix QualiPlate Kit for Cry34Ab1 is designed for the qualitative laboratory detection of Cry34Ab1 protein in corn leaf or seed. It may be used to detect the presence of the protein in corn products including HERCULEX™ RW and HERCULEX XTRA.

How the Test Works

This Kit is a “sandwich” Enzyme-Linked ImmunoSorbent Assay (ELISA). In the test, corn sample extracts are added to test wells coated with antibodies raised against Cry34Ab1 protein. Any Cry34Ab1 present in the sample extract binds to the antibodies and is then detected by addition of enzyme (horseradish peroxidase)-labeled Cry34Ab1 antibody.

After a simple wash step, the results of the assay are visualized with a color development step. Color development increases with increasing Cry34Ab1 sample concentration.

Lighter color = Low concentration
Darker color = High concentration

How the Kit Performs

This Kit is a strictly qualitative (yes/no) assay. Samples are interpreted in comparison to a Positive Control provided in each kit. Instructions for interpreting results based upon this control start on page 3. It is recommended that the user also prepare known negative and positive seed or leaf samples be run in every assay as controls, in addition to the kit Positive Control.

Preparation of Solutions

Note: To handle bulk packaged Cry34Ab1 Enzyme Conjugate and Substrate, pour off 6 milliliters of Conjugate and 11.5 mL Substrate per plate to be run each day. Use a multiple-channel pipette to dispense. Do not pour excess Substrate back into the reagent bottle.

Wash Buffer (PBS-0.05% Tween 20, pH 7.4): May be purchased in 1L dry packets from Sigma Chemicals, Cat#P-3563, or prepared from salts on site. Store at controlled ambient temperature for up to one week, then discard. At least 100 mL/plate is needed, plus additional for preparation of Extraction Buffer.

Extraction Buffer: Add 0.5 mL Tween 20 to 100 mL Wash Buffer. Store refrigerated when not in use; warm to room temperature prior to assay.

Stop Solution (1 N Hydrochloric acid (HCl): Prepare by adding 83 mL of concentrated HCl (36.5-38.0%) to 917 mL of distilled or deionized water; work in a fume hood and use proper protective gear. This reagent may be stored at room temperature for 2 years.

Sample Preparation

Single Corn Seed Samples:


1. Crush seeds: Seeds may be crushed by any number of methods, from hammers or pliers in a bag or tube, to 48-well seed crushers, to bead-beater type grinders.



Prepare wash buffer and grain extraction solutions



Allow all reagents to reach room temperature before beginning



Add controls and samples

Whatever the method used, take extreme care not to cross-contaminate between seed samples.

CAUTION: Cry34Ab1 protein is expressed at high concentrations in corn seed, so there is serious potential for cross-contamination between samples during seed crushing. Use the utmost care to avoid this. Cleaning the cutting/crushing surfaces with an alcohol-soaked pad between samples is recommended.

2. Add 1 mL of Extraction Buffer to each crushed seed. Mix for at least 30 seconds. For best results, allow to extract for an hour, mixing again at the end of that time. If seeds are thoroughly crushed, this extraction time can be reduced. Allow extracts to settle completely. Dispensing particles into the test plate can cause false positive results.

Single Corn Leaf Punch Samples:

1. Take a single leaf punch of approximately 5 to 10 millimeters diameter, using a paper punch or micro-tube cap. Mash the leaf tissue with a pestle matched to the micro-tube, or disrupt via another method. The extraction efficiency of whatever method used will vary proportionately with the amount of tissue disruption performed.
2. Add 0.25 to 0.5 mL of Extraction Buffer per leaf punch. Mash again with pestle or mix vigorously for at least 30 seconds, and allow solids to settle. Dispensing particles into the test plate can cause false positive results. Take extreme care not to cross-contaminate between leaf samples. Well-extracted leaves should result in a green cloudy extract. Note the presence of any extracts that appear clear and/or colorless – these may not have extracted properly and assay data would be invalid. For best results, extract another representative sample.

How to Run the Assay


- Read all of these instructions before running the kit.
- Allow all reagents to reach room temperature before beginning (at least 30 minutes with un-boxed plates and reagents at room temperature - do not remove strips from bag with desiccant until it has warmed up).
- Organize all Control and sample extracts, and pipettes so that Step 1 can be performed in 15 minutes or less. The use of a multi-channel pipette is strongly recommended for all reagent additions.
- Use the well identification markings on the plate frame as a guide when adding the samples and reagents. For this qualitative assay, duplicate wells of the Extraction Buffer Blank (BL) and the Positive Control (PC), along with 92 sample extracts and user-supplied controls (S) in single wells may be run on one plate. (See the Qualitative Assay Example Plate Layout - Figure 1a).

Procedure

1. Add **50 μ L** of **Extraction Buffer Blank (BL)**, **50 μ L** of **Positive Control (PC)**, and **50 μ L** of each **sample and user-prepared control extract (S)** to their respective wells, as shown in the Example Plate Layout (Figure 1A).

NOTE: In order to minimize setup time it is strongly recommended that a multi-channel pipette be used in steps 1, 5, 8, and 10.

2. Thoroughly mix the contents of the wells by moving the plate in a rapid circular motion on the benchtop for a full 20-30 seconds. Be careful not to spill the contents!
3. Cover the wells with tape or Parafilm to prevent evaporation and incubate at ambient temperature for **15 minutes**. If an orbital plate shaker is available, shake plate at 200 rpm.



Mix plate



Incubate



Bottle Wash method



Strip Plate Wash option

4. After incubation, carefully remove the covering and vigorously shake the contents of the wells into a sink or other suitable container. Flood the wells completely with **Wash Buffer**, then shake to empty. Repeat this wash step three times. Alternatively, perform these four washes (300 μ L/well) with a microtiter plate or strip washer. Slap the inverted plate on a paper towel to remove as much liquid as possible.
5. Add **50 μ L Cry34Ab1 Antibody-Enzyme Conjugate** to each well. Thoroughly mix the contents of the wells, as in step 2.
6. Cover the wells with tape or Parafilm to prevent evaporation and incubate at ambient temperature for **45 minutes**. If an orbital plate shaker is available, shake plate at 200 rpm.
7. Wash the plate as described in step 4.
8. Add **100 μ L of Substrate** to each well.
9. Thoroughly mix the contents of the wells, as in step 2. Cover the wells with new tape or Parafilm and **incubate for 15 minutes at ambient temperature**. Use orbital shaker if available.

Caution: Stop Solution is 1.0N Hydrochloric acid. Handle carefully.

10. Add **100 μ L of Stop Solution** to each well and mix thoroughly. This will turn the well contents yellow.

NOTE: Read the plate within 30 minutes of the addition of Stop Solution.

How to Interpret the Results

Spectrophotometric Measurement

1. Set the wavelength of the microtiter plate reader to 450 nanometers (nm). (If it has dual wavelength capability, use 600, 630 or 650 nm as the reference wavelength.)
2. Set the plate reader to blank on the Extraction Buffer Blank wells (this should automatically subtract the mean optical density (OD) of the Blank wells from each control and sample OD). If the reader cannot do this, it must be done manually.

General test criteria:

The mean OD of the BLANK wells should not exceed 0.15.

The mean, blank-subtracted OD of the Positive Control wells should be at least 0.2.

The coefficient of variance (%CV) between the duplicate Positive Control wells should not exceed 15%:

$$\%CV = \frac{\text{std. deviation of OD's}}{\text{mean Pos.Ct. OD}} \times 100$$

If the results of an assay fail to meet these criteria, consult EnviroLogix' Technical Service for suggestions on improving the test when you repeat the assay.

Calculate the Positive Control Ratio

Divide the OD of each sample extract by the mean OD of the Positive Control wells. This number is the "Positive Control Ratio".



*Read plates in a Plate Reader
within 30 minutes of the
addition of Stop Solution*

Interpret the Qualitative Results

Single Corn Leaf and Seed samples:

If the Positive Control Ratio calculated for a sample is less than 0.5, the sample does not contain Cry34Ab1 protein.

If the Positive Control Ratio of a sample is greater than or equal to 0.5, the sample does contain Cry34Ab1 protein.

Leaf and seed samples are by their nature either 100% positive or 100% negative. Any low level positive results from single seed or leaf samples must be due to either some form of sample cross-contamination (stray particles or dust from corn, leaf residue on leaf punch, etc.) or can be caused by transfer of particulate matter from leaf or seed extracts into the assay wells. If there is any question of the latter occurring, re-extraction and re-testing is recommended.

Figure 1a. Example of a typical Qualitative assay setup.

	1	2	3	4	5	6	7	8	9	10	11	12
A	BL	S7	S15	S23	S31	S39	S47	S55	S63	S71	S79	S87
B	PC	S8	S16	S24	S32	S40	S48	S56	S64	S72	S80	S88
C	S1	S90	S17	S25	S33	S41	S49	S57	S65	S73	S81	S89
D	S2	S10	S18	S26	S34	S42	S50	S58	S66	S74	S82	S90
E	S3	S11	S19	S27	S35	S43	S51	S59	S67	S75	S83	S91
F	S4	S12	S20	S28	S36	S44	S52	S60	S68	S76	S84	S92
G	S5	S13	S21	S29	S37	S45	S53	S61	S69	S77	S85	BL
H	S6	S14	S22	S30	S38	S46	S54	S62	S70	S78	S86	PC

Precautions and Notes

- Store all QualiPlate Kit components at 4°C to 8°C (39°F to 46°F) when not in use.
- Do not expose QualiPlate Kit components to temperatures greater than 37°C (99°F) or less than 2°C (36°F).
- Allow all reagents to reach ambient temperature (18°C to 27°C or 64°F to 81°F) before use.
- Do not use kit components after the expiration date.
- Do not use reagents or plates from one QualiPlate Kit with reagents or plates from a different QualiPlate Kit.
- **Do not expose Substrate to sunlight** during pipetting or while incubating in the test wells.
- The assay has been optimized to be used with the protocol provided in the kit. Deviation from this protocol may invalidate the results of the test.
- As with all tests, it is recommended that results be confirmed by an alternate method when necessary.
- Observe any applicable regulations when disposing of samples and kit reagents.
- Use caution to prevent sample-to-sample cross-contamination with samples, fluids, or disposables.



**For Technical Support
Contact Us At:**

EnviroLogix
500 Riverside Industrial
Parkway
Portland, ME 04103-1486
USA

Tel: (207) 797-0300
Toll Free: 866-408-4597
Fax: (207) 797-7533

e-mail:
info@envirologix.com

website:
www.envirologix.com



LIMITED WARRANTY

EnviroLogix Inc. ("EnviroLogix") warrants the products sold hereunder ("the Products") against defects in materials and workmanship when used in accordance with the applicable instructions a period not to extend beyond a product's printed expiration date. If the Products do not conform to this Limited Warranty and the customer notifies EnviroLogix in writing of such defects during the warranty period, including an offer by the customer to return the Products to EnviroLogix for evaluation, EnviroLogix will repair or replace, at its option, any product or part thereof that proves defective in materials or workmanship within the warranty period.

ENVIROLOGIX MAKES NO OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. The warranty provided herein and the data, specifications and descriptions of EnviroLogix products appearing in EnviroLogix published catalogues and product literature are EnviroLogix' sole representations concerning the Products and warranty. No other statements or representations, written or oral, by EnviroLogix' employees, agents or representatives, except written statements signed by a duly authorized officer of EnviroLogix Inc., are authorized; they should not be relied upon by the customer and are not a part of the contract of sale or of this warranty.

EnviroLogix does not warrant against damages or defects arising in shipping or handling, or out of accident or improper or abnormal use of the Products; against defects in products or components not manufactured by EnviroLogix, or against damages resulting from such non-EnviroLogix made products or components. EnviroLogix passes on to customer the warranty it received (if any) from the maker thereof of such non-EnviroLogix made products or components. This warranty also does not apply to Products to which changes or modifications have been made or attempted by persons other than pursuant to written authorization by EnviroLogix.

THIS WARRANTY IS EXCLUSIVE. The sole and exclusive obligation of EnviroLogix shall be to repair or replace the defective Products in the manner and for the period provided above. EnviroLogix shall not have any other obligation with respect to the Products or any part thereof, whether based on contract, tort, strict liability or otherwise. Under no circumstances, whether based on this Limited Warranty or otherwise, shall EnviroLogix be liable for incidental, special, or consequential damages.

This Limited Warranty states the entire obligation of EnviroLogix with respect to the Products. If any part of this Limited Warranty is determined to be void or illegal, the remainder shall remain in full force and effect.

This test kit has been validated and approved by Dow AgroSciences for detection of the Cry34Ab1 protein expressed in corn products including HERCULEX RW and HERCULEX XTRA.

Parafilm is a registered trademark of American Can Corporation
HERCULEX is a trademark of Dow AgroSciences LLC
Tween is a registered trademark of Uniqema, a business unit of ICI Americas Inc.
EnviroLogix, the EnviroLogix logo, and QualiPlate are trademarks of EnviroLogix Inc.

© EnviroLogix 2007